

Programmable Low Temperature Precision Bath (Pro-cool bath)

NCB-3100 · 3200 · 3300




NCB-3100

NCB-3200

NCB-3300

Operation unit



Many program functions for various purposes

- Wide range temperature control with the precision at $\pm 0.03^{\circ}\text{C}$
- The temperature trend data can be logged and analyzed feeding out to PC.
- Provide programs of up-rising and down falling with precise control. Seven program patterns (Max. 10 segments per pattern) are available, and patterns can be spliced one another.
- The post operation at the end of the program, can be selected from three endings, Terminate after the end, preservation of the end temperature and back to start.

Model	NCB-3100		NCB-3200	NCB-3300	
Product code No.	112450		112460	112470	
Performance	Temperature control range / precision		$-30^{\circ}\text{C} \sim 80^{\circ}\text{C} / \pm 0.03^{\circ}\text{C}$		
	Cooling capacity (at liquid temperature)	20°C	370W (318kcal/h)	480W (412kcal/h)	610W (524kcal/h)
		-10°C	220W (189kcal/h)	350W (301kcal/h)	420W (361kcal/h)
	Temperature fall slope	1.1°C or below/min. (40→0°C)		0.3°C or below/min.	0.3°C or below/min.
0.5°C or below/min. (5→-15°C)					
0.4°C or below/min. (-15→-30°C)					
External circulation ability	Maximum discharge 6.7L/min (50Hz, water), Maximum pump head 2.9m (50Hz, water)				
Function	Temperature control (sensor)	P.I.D control, SSR output (platinum resistance temperature detector JPt 100Q)			
	Temperature setting / Display	Numeric keypad input type, Least input digit: 0.01°C, Digital display			
	Programming functions	1. Fixed value operation (Normal mode), 2. Auto start, 3. Auto stop (1 min. to 99 hrs. and 59min.), 4. 7 patterns of user programs (Max 10 segments / pattern, 1 min. to 99 hrs. and 59 min / segment)			
	Safety functions	Self-diagnostic function (water level decrease, sensor malfunction, SSR defect, overheating, heater breaking, cooling defect, power failure alarm, program slope defect), Breaker for electrical leakage and overcurrent, Refrigerator overload relay			
External input / output	RS-232C connecting terminal, External temperature sensor terminal				
System	Heater	1.2kW (SUS 316L)			
	Refrigerator/Coolant	Air cooling type, Output 450W / R404A	Air cooling type, Output 400W / R407C	Air cooling type, Output 400W / R404A	
	Cooling coil	SUS304			
	Inner bath stirring	Water jet type			
Specification	Nozzle for external circulation	Closed loop circulation, Discharge port: OD 10.5 mm hose nozzle (Rc1/4, with stop valve), Return port: OD 10.5 mm house nozzle (Rc1/4)			
	Material of bath	SUS304			
	Bath inside dimensions(mm) / capacity	250W×295D×180H / 13L	320W×405D×180H / 23L	430W×517D×180H / 41L	
	Bath effective dimensions(mm)	215W×135D×160H	280W×240D×160H	400W×354D×160H	
Surrounding temperature range	5-30°C				
Outside dimensions	620(635)W×390(440)D×464H	393W×511D×793H	510(525)W×627(667)D×898H		
Weight	about 46 kg	about 48 kg	about 59 kg		
Power input / voltage	20A, 2kVA / AC 100V, 50/60Hz with no plug				

※The above performance is yielded at RT 20°C, at liquid temperature of 40°C, at rated power supply voltage of 50 Hz, using water, at no load. ※Thermal media specialized for low temperatures range is required in case of the usage below 10 °C. The temperature control range is -30 to 70°C at surrounding temperature of 20°C without putting a lid on. (NCB-3300)

※Dimensions shown in parentheses include projections.

Power plug is not supplied.

POINT Interface with PC & External Circulation




Self-diagnosis function is on board.
Stop the control and inform by alarm upon detection of abnormality in the unit. Stop valve nozzle is equipped for extra-tank circulation. Temperature control of external units is also possible by using optional temperature sensor.

External Temperature Sensor STP-200
Product code No.113620

Operation pattern

■ User program



As User programs, we have A type (3 patterns) which stops the control at the end of the program and B type (4 patterns) which keeps the final temperature even after ending the program by stand-by function.

For repeated experiments, if you use B type, you can keep the temperature by driving SKIP/END function during stand-by operation and changing back to the default. The program runs repeatedly when you set a new sample and press the key again.

Options


Clamp Set:

Pole, holder and clamp are included so that you can fix flask etc.
Product code No.164500

POINT Remote Control & Data Processing by Using Supplied Exclusive Software

Simplified measuring/controlling software ESMon2018

You can control temperature remotely through RS-232C communication cable using PC and process the trend data.
※Only normal mode is available, and the program mode is not supported.




ESMon2018 Operation Environment

Supported OS: Microsoft Windows®7/10
CPU : Pentium or above (recommend Pentium III 800MHz or above)
Hard disk : 15MB or more of free space is required.
More space is necessary when you record the data.
Memory : 128MB or more is recommended though it depends on OS.
Communication Interface: Serial interface (COM1 etc.) is necessary. Please use
Please use commercial "USB-serial conversion cable" etc. when you use PC without serial interface.


Data

■ Cooling Capacity Curve




■ Heating & Cooling Time

NCB-3200



■ Heating & Cooling Time

NCB-3100



■ Heating & Cooling Time

NCB-3300

