


Spray Dryer

Model SD-1010



Model SD-1010

Control panel



External output



	Output range	Electrical signal
Inlet temperature	50 to 250°C	0 to 10mV
Outlet temperature	50 to 250°C	0 to 10mV
Dry air volume	0 to 1m ³ /min	0 to 10mV

Optional set

Model SD-1010 with recorder (Cat. No. 275190)

Recorder (μR-10000) with adjusted input range for dry air volume and temperature at 2 points is included.

For polymers, pharmaceuticals, organic and inorganic compounds, and new materials

- Liquid samples are directly converted into fine powder. There is no need for filtration or concentration as a pre-process, or pulverization as a post-process.
- Instantaneous drying requires a short processing time. No high-temperature load is applied to the sample.
- Experiments can be easily performed at the laboratory level since the sample can be dried even with small amounts of solids on the order of grams.
- The unit is equipped with a digital display of the inlet and outlet temperatures and the amount of drying air, as well as a recorder output for setting, checking, and recording experimental data.
- The nozzle part of the spray air line is automatically cleaned by time setting. Clogged nozzles in the sample liquid line can be cleaned by pressing the needle. A cleaning port is provided on the evaporation tube, allowing easy removal of crystals adhering to the tip of the nozzle.
- In order to prevent liquid accumulation and sample sedimentation, the layout is designed to take into account the position of the nozzle for spraying, the sample stirrer and the pump.
- To prevent disconnection or damage to the liquid delivery tube due to clogging of the spray nozzle, a discharge pressure monitoring function is provided for the liquid delivery pump.
- Equipped as standard with a magnetic stirrer, inlet (suction) filter, and outlet (air) filter to prevent sedimentation and separation of samples.
- Easy disassembly, cleaning, and assembly of equipment.
- Data (operating conditions and recovery rates from various samples) are available on our website.

Model	SD-1010
Cat. No.	275178
Type	Spray nozzle (Two-fluid nozzle)
Moisture evaporation	Max. 1500mL/h
Temp. control range/accuracy	40 to 200°C (at inlet port) / ±1°C
Heated air volume control range	0.2 to 0.75m ³ /min
Spray air pressure control range	20 to 250kPa (0.2 to 2.5kg/cm ²)
Pump flow rate range	150 to 1700mL/h
Stirring rotation speed/capacity	100 to 1000rpm/50mL to 2L (water)
Safety functions	Evaporation tube protection cover, earth leakage, overcurrent breaker, fixed temp. overheat protector, automatic recovery from power failure, self-diagnostic function (temp. abnormality, air volume abnormality, heater disconnection, abnormal dispensing pressure)
Spray air cleaning function	Automatic operation by time interval
Additional functions	Spray air line valve ON/OFF switch
External output	0 to 10mV (inlet/outlet port temp., dry air volume)
Heater	3kW
Liquid dispensing pump	Peristaltic pump
Spray nozzle	Two-fluid nozzle (nozzle hole ID 0.71mm)
Spray air pump	Air compressor (option)
Evaporation tube	Hard glass
Cyclone	Hard glass
Powder collection vessel	Hard glass (capacity 600mL)
Sample stirrer	Brushless DC motor 6W, cobalt magnet
Suction blower	Commutator motor, Max. 0.75m ³ /min variable flow type
Tube diameter	ID3.15 x OD5.2mm
Spray air connecting port diameter	ID4 x OD6mm soft urethan tube union
Spray air pressure	Pressure 294kPa (3kg/cm ²), flow rate 25L/min or more
Exhaust connecting port diameter	Exhaust port OD50mm
Interval timer	OFF or 1 to 20min (automatic intermitted time)
Ambient temperature range	5 to 35°C
Outer dimension (mm)/Weight	700W x 620D x 1500H / approx. 110kg
Electrical connection	19.1A, 4.2kVA / AC220V single phase 50/60Hz

※Specifications are based on room temperature at 20°C, rated power supply voltage, 50 Hz, tap water, and use of silicone tubing.

✳ Power plug is not included.


Accessories	
	Silicone tube ID3.15 x OD5.2mm x 2m Exhaust hose 2.5m, recorder output cable x 3pcs Air piping tube ID4mm x OD6mm x 5m, Air piping connector (1/4B)



Designed for safety and ease of use

Point

Safety design



Spray dryer, Model SD-1010 is equipped with abundant safety functions, including self-diagnostic alarms for abnormal temperature, airflow, heater disconnection, and pumping pressure, leakage/over-current breaker, overheat protector and manual recovery in the event of power failure. In addition, a rubber tray is provided at the tube piping for safety design.

Streamlined layout of tubing lines




The tubing line is laid out in a streamline position and at the shortest possible distance, preventing liquid accumulation and sample sedimentation.

Abnormal pressure detection function




In order to prevent the dispensing tube from being disconnected or damaged due to a clogged nozzle, the pressure sensor is set in the dispensing and pressure lines to issue a pressure abnormality alarm.



Safety design

The device is equipped with a protective cover to prevent accidental touch on the evaporation tube which becomes hot during operation.

Point Easy maintenance



Equipped with a stand to prevent evaporation tubes from falling

Simple disassembly, cleaning and assembly

The conventional troublesome disassembly, cleaning, and reassembly can be performed smoothly. The cyclone can be disassembled, and position alignment marks are engraved for easy setting of the evaporation tube. Filters on the inlet and outlet sides can be easily replaced through the maintenance doors on the left and right sides of the main unit.

Suction filter (inlet)



Prevents external contamination.

Air filter (outlet)




Protects the suction blower.




Flow sheet

Model SD-1010



Caps for nozzles



Caps (Hole diameter mm)		Nominal diameter	Spray air Flow rate at P pressure [kPa (kg/cm ²)] (L/min)	Sample type	Cat. No.
For liquid	PF-1650-SS (0.41)	1A	68.67 (0.7) → 1.3 to 147.15 (1.5) → 17.0	Homogeneous solution	120710
For air	PA-64-SS (1.6)				120740
For liquid	PF-2050-SS (0.51)	1	68.67 (0.7) → 13.3 to 147.15 (1.5) → 20.0	Slurry concentration Approx. 5%	120720
For air	PA-64-SS (1.6)				120740
For liquid	PF-2050-SS (0.51)	2A	68.67 (0.7) → 23.0 to 147.15 (1.5) → 36.0	Slurry concentration Approx. 5%	120720
For air	PA-70-SS (1.75)				120750
For liquid	PF-2850-SS (0.71)	2 (standard)	68.67 (0.7) → 19.3 to 147.15 (1.5) → 31.0	Slurry concentration Approx. 5 to 20%	120730
For air	PF-70-SS (1.75)				120750

Needles (for spray air nozzle)

Hole diameter	For 0.41mm	For 0.51mm	For 0.71mm
Cat. No.	180250	180260	180270

Exhaust hose

2.5m Cat. No. 179980
*Standard accessories

10m Cat. No. 179990



Silicone tube

(ID3.15 x OD5.2mm x 5m)
Cat. No. 125510

Hot plate unit
Model HTP-1000

Ultimate temperature
Max. 180°C (plate surface)
Stirring table dimensions
φ 135mm
Cat. No. 275300
External temperature sensor
Cat. No. 265680

Vessel set with drain port
Cat. No. 275640


Set configuration

① Vessel with drain port
Cat. No. 275210

② Two-way cock set
Cat. No. 275220

③ Pole for SD-1010
Cat. No. 275200

④ Holder for vessel with drain
port
Cat. No. 275230



Suction filter (inlet)
Cat. No. 179950
*Standard accessories



Air filter (outlet)
Cat. No. 179960
*Standard accessories



Separator / powder collector (600mL)
Cat. No. 146750
*Standard accessories
Powder collector (1.2L)
Cat. No. 147940

*Resin-coated EYELACOAT® bottles that can
reduce shattering of glass are also available.
Please contact us for more details.

Recorder



Product name	Recorder
Model	μR-10000 (for SD-1010)
Cat. No.	219260
Measurement points	6 points/6 ribbon colors
Measurement cycle	6 points/10sec (fastest)
Input	DCV: 5 types, TC: 1 type, RTD: 2 types

Air compressor



Product name	Air compressor	
Model	0.2LE-8SBA	0.4LE-8SB
Discharge air volume (50/60Hz)	20/24L/min	42/49L/min
Tank capacity	30L	
Control pressure (kgf/cm ²)	3.5 to 5.0	5.0 to 7.5

*Please source an equivalent model as a base locally.