


Temperature Ramp Incubator

**5 Chambers with Individual Temperature Control
LED Illumination for Day/Night Cycle Operation**

Model MTI-1000/1100



Precise individual temperature control

The Temperature Ramp Incubator has the function of five (5) independent temperature-controlled chambers, allowing users to set up different environment conditions in each chamber.

[Example].....

The same seeds were exposed to five different temperatures, simulating 14-hour day and 10-hour night, and their germination was observed.

• Conditions

Illumination: 4 lamps (15,000lx), 14 hours ON, 10 hours OFF with 24-hour cycle timer.

Chamber	Temperature
5	32°C
4	28°C
3	24°C
2	20°C
1	16°C




Equipment used
MTI-1100

Standard 40W x 4 LED lamps (cool white)

The illumination intensity can be adjusted by turning on 1 to 4 LED lamps, and a 24-hour cycle timer allows for creating day and night environments. Cool white LED lamps provide brightness of up to 15,000lx and a long life of approximately 40,000 hours.


Relative spectral radiant intensity of cool white LED lamps



Equipped with 2-step program

The model MTI-1100 allows you to set 2-step program to reproduce day and night temperatures. Day/night tests can be performed by switching the cool white LED lamps on and off using the ON-OFF timer.

2-step programmable temperature control



The temperature can be controlled in two steps to turn the LED lamp on and off. Day/night operation is possible based on temperature and illuminance.

Low GWP Refrigerant R448A (GWP value 1386)

※GWP stands for Global Warming Potential. It represents the capacity of a greenhouse gas to contribute to global warming compared to carbon dioxide, which has a GWP value of 1.



Chamber illumination

Cool white straight LED lamps
 1 lamp ON: Approx. 4,000lx
 2 lamps ON: Approx. 8,000lx
 3 lamps ON: Approx. 12,000lx
 4 lamps ON: Approx. 15,000lx




Power supply
 100V 15A

It can be powered from a household outlet.

Equipped with standard casters and adjusters

Simple measurement and control software (ESMon)



Temperature data for each chamber can be obtained on PC, and the data can be recorded and processed.

In addition, the model MTI-1100 can output the temperature of each chamber as analog (0 to 5 VDC), so it can be connected to a data logger or recorder.

Specifications

Model	MTI-1000	MTI-1100
Cat. No.	283330	283340
Temp. control range	5 to 50°C	
Temp. control accuracy	Chamber 1: ±2.5°C~, Chamber 2: ±1.5°C~, Chamber 3 to 5: ±1.0°C~	
Number of chambers	5 chambers	
Illumination	Cool white straight tube LED lamp 40W x 4 pcs	
Temp. control method	Constant temperature control	2-step temperature control
Temp. control	Chamber 1 (bottom level) Refrigerator ON-OFF control (NORM. mode), Refrigerator ON-OFF + timer control (INT. mode), Other chamber heater control (P.I.D.)	
Illumination time	24-hour timer (minimum setting 15 min.)	
Heater	600W	
Refrigerator/refrigerant	Air-cooled, 400W output, R448A (GWP value 1386)	
External input/output	USB communication terminal, simplified software ESMon for model MTI (accessory)	
Shelf load capacity / Material	Equal distribution Max.25kg/shelf, Tempered glass	
Inner outlet	1 port in each chamber (5 total 0.5A max.)	
Inner chamber dimensions / Capacity	400(375)W×470D×195H (mm)/Approx. 37L	
Ambient temperature range	5 to 30°C	
External dimensions / Weight	560W x 670(711)D x 1780H(mm)/Approx. 160kg	
Power input/Rated power supply	15A, 1.5kVA 100VAC 50/60Hz	

※Performance is measured in NORM. mode, room temperature 20°C, rated power supply voltage 50Hz, no illumination, and no load.

※Illuminance is measured from the center of the chamber facing rearward. Illuminance decreases as LED lamps deteriorate.

※Power input includes the capacity of the outlet in the chamber (0.5A).


※Temperature setting of each chamber is subject to conditions.

A shaker can be installed for shaking culture.

One unit of shaker can be installed in the chamber 1. The compatible shaker models are MMS-120H, 220H, 320, and 420.

※The shaking speed of the shaker must be within 150 rpm for rotary shaking and 100 rpm for reciprocal and figure-8 shaking.

※The distance between the top surface of the shaker and the shelf board above is 75mm.



Options

Recorder

Model μR-10000

Cat. No. 217060

※Only for MTI-1100

Recorder output cord (1 pc)


Cat No. 147570

※Five (5) output cords are required to record temperatures in the chamber 1 to 5.


Data

Temperature control accuracy

Conditions: Room temperature 20°C, 100VAC 50Hz, LED off, no load in the chamber, temperature settings 9, 18, 27, 36, 45°C, sensor position in the center of each chamber



Temperature setting range



This graph shows the possible settings for the chamber 5 and between each chamber, assuming the same temperature difference across the five chambers. When the room temperature is 5°C, control of the orange section becomes unavailable. Similarly, when the room temperature is 30°C, control of the purple section becomes unavailable. A temperature difference must be maintained between each chamber. (Example) If the chamber 5 is set to 30°C, the temperature difference between each chamber can range from 2°C to 6°C.

TOKYO RIKAKIKAI CO., LTD.

EYELA

<https://eyelaworld.com>

TN Koishikawa Bldg.
 1-15-17 Koishikawa
 Bunkyo-ku, Tokyo
 112-0002 Japan

Tel: +81-3-6757-3378
 Fax: +81-3-3868-6571
 E-mail: contact@eyelaworld.com

※The appearance and specifications of the products are subject to change without notice.

2024.06_K5