

Compact Water Bath NTT-2000

Mini bath/shaker for various uses. Economical and low power consumption.



NTT-2000

Features

- Compact (200Wx325Dmm) and low power consumption (4A).
- Mild stirring realizes low water level incubation as low as 40mm water surface level.
- Timer function is provided (Auto start, Auto stop).
- Safety functions include self diagnosis, float switch and independent overheat protector.

Specifications

Product Name	Water Bath
Model	NTT-2000
Cat. No.	195400/195408
Temp. control range	RT + 5 ~ 70°C
Temp. control accuracy	±0.07°C
Temp. control	Microprocessor P.I.D. control, SSR output
Safety functions	Self diagnosis (Temp. limit, Heater brake, Sensor brake), Empty heating protector, Independent overheat protector, Fuse
Additional functions	Alarm buzzer, Auto tuning
Heater	300W (SUS 316L Water only)
Dimensions (mm) & Weight	200W x 105D x 207H 2.5kg
Power source	400VA AC100V/880VA AC220V

Option



Bath cover SBK-2
Cat. No. 196200

Shaker SS-1000

Exclusive shaker for use with NTT-2000

By placing NTT-2000 on top of the SS-1000, the system functions as small stirring water bath.



SS-1000

NTT-2000 + SS-1000

Specifications

Product Name	Shaker
Model	SS-1000
Cat. No.	195460
Shaking method	Reciprocating shaker
Shaking speed	0 ~ 200 reciprocate/min
Shaking width	25mm
Dimensions (mm) & Weight	200W x 469D x 135H 6.6kg
Power source	40VA AC100V

Rotary Evaporator 1

Vacuum Pump & Aspirator 2

Organic Synthesizer 3

Freeze Dryer - Strip Dryer 4

Shaker 5

Stirrer 6

Tube Pump 7

Low Temperature Oscillator 8

Cold Trap 9

Constant Temperature Bath 10

Incubator 11

Oven 12

Fermentor 13