




Rotary evaporator combination guide

You can freely select the configuration you want according to your needs




Type E condenser
Vertical double spiral

Adapter integrated type, suitable for installation inside a fume hood




Type S condenser
Diagonal double spiral

Highly efficient type for low to high boiling point substances



Type V condenser
Vertical double spiral + adapter


Accommodates a wide range of flask sizes, from small to large volumes




Type T condenser
Dewar flask type

Suitable for low-boiling-point substances using ice or dry ice-acetone

► Driving unit



Standard type
Manual jack
Model N-1300
(Type E · V · S)




High-performance type
Motorized jack
Model N-1210B
(Type V · S · T)

► Bath



For water
SB-1300






For water·oil
OSB-2200

► Rotary evaporator N-1300 series

Condenser		Type E condenser		Type V condenser		Type S condenser
Model		E-29	EF-29 EYELA COAT®	V-29	VF-29 EYELA COAT®	S-29
Cat. No.		266310	—	266290	—	266270
		N-1300E		N-1300V		N-1300S
Driving unit		N-1300E	N-1300EF	N-1300V	N-1300VF	N-1300S
Model	N-1300					
Water bath		N-1300E-W	N-1300EF-W	N-1300V-W	N-1300VF-W	N-1300S-W
Model	SB-1300					
Water · oil bath		N-1300E-WB	N-1300EF-WB	N-1300V-WB	N-1300VF-WB	N-1300S-WB
Model	OSB-2200					

► Rotary evaporator N-1210B series

Condenser		Type V condenser		Type S condenser	Type T condenser
	Model	V-29	VF-29 EYELA COAT®	S-29	T-29
	Cat. No.	269990	—	266270	270030
					
					
Driving unit		N-1210BV	N-1210BVF	N-1210BS	N-1210BT
Model	N-1210B				
					
Water bath		N-1210BV-W	N-1210BVF-W	N-1210BS-W	N-1210BT-W
Model	SB-1300				
					
Water · oil bath		N-1210BV-WB	N-1210BVF-WB	N-1210BS-WB	N-1210BT-WB
Model	OSB-2200				

► N-1300 · N-1210B classification

Jack type		Manual jack			Motorized jack		
Bath type		Without bath	Water bath	Water · oil bath	Without bath	Water bath	Water · oil bath
Condenser Type E	Vertical double spiral	N-1300E	N-1300E-W	N-1300E-WB			
Type EF	Vertical double spiral EYELA COAT®	N-1300EF	N-1300EF-W	N-1300EF-WB			
Type V	Vertical double spiral + adapter	N-1300V	N-1300V-W	N-1300V-WB			
Type VF	Vertical double spiral + adapter EYELA COAT®	N-1300VF	N-1300VF-W	N-1300VF-WB	N-1210BVF	N-1210BVF-W	N-1210BVF-WB
Type S	Diagonal double spiral	N-1300S	N-1300S-W	N-1300S-WB	N-1210BS	N-1210BS-W	N-1210BS-WB
Type T	Dewar flask type				N-1210BT	N-1210BT-W	N-1210BT-WB